


EyeZoom™

The first and only adjustable magnification loupe, EyeZoom offers three magnification levels ranging from 3x-5x.

Developed to assist healthcare professionals view an operating site from multiple perspectives, the loupe conforms to fit each user's individual preference.

Co-engineered with Konica Minolta®, the patent-pending technology provides edge-to-edge clarity and high-definition resolution. With an ideal symmetry between style and performance, the lightweight magnesium bezel leaves a lightweight feel and an aesthetically pleasing design

Konica Minolta® is a registered trademark of Konica Minolta, Inc. Used with permission.



Mounting Options:	Through-the-Lens
Magnification Power:	3x, 4x and 5x
Field Width:	4" (10.2 cm) at 3x, 3.25" (8.3 cm) at 4x, 2.5" (6.35 cm) at 5x
Field Depth:	4" (10.2 cm) at 3x, 3.25" (8.3 cm) at 4x, 2.5" (6.35 cm) at 5x
Working Distance:	Customized to each user
Declination Angle:	Customized to each user
Weight*:	3.4 oz (97.6 g)
Loupe package:	Side shields Head strap Personalized storage case Optical screwdriver Cleaning cloth
Frame Options:	Legend™ Rave™ Rydon™ Victory™ XV1™

All optical systems are made with Grade A fine annealed glass lenses fused with anti-scratch and anti-reflective coatings. All carrier lenses are made with High Index polycarbonate material.

* Measured on Orascoptic Legend frame

BERNER
Medical

Bröderna Berner Handels AB, Division Berner Medical, Box 50132, 202 11 Malmö
tfn 040-680 68 00, info@bernermedical.se, www.bernermedical.se