

# Endeavour™ XL

The Endeavour XL headlight provides a precise beam of light that helps illuminate critical details that can often be overlooked by the naked eye.

With 25% longer run-time than the standard Endeavour headlight, the Endeavour XL provides up to 24 hours of continuous battery life.

The addition of multiple light intensities up to 85 lumens and 25% brighter light output than the Endeavour, the headlight allows the user to select the optimal light output for each procedure.

Constant-current technology also ensures the headlight's light output will not diminish during the run-time of the battery pack.

The Endeavour XL headlight comes with a three year warranty on both the cable and LED and one year warranty on the battery pack.



## Headlight Specifications

Weight:	0.38 oz (10.8 g)
Light Intensity (Low/Medium/High):	53/68/85 lumens
Intensity Control:	Capacitive Touch up to 85 lumens
Color Temperature:	6500° K
Focus:	2.5" (6.35 cm) spotlight at 12" (30.5 cm) focused beam
Composite Filter:	Curing filter included

## Battery Specifications

Dimension:	3.6" (9.2 cm) x 2.1" (5.5 cm) x 1.1" (2.8 cm)
Est. Battery Life	10 to 24 hours depending on light setting
Est. Charging Time:	up to 8 hours (can be used while recharging)
Type:	Lithium Ion

