

Pulmonary[®]

BiTrac Select™

ONE Mask. MULTIPLE uses.

Pulmonary's BiTrac Select series of noninvasive ventilation masks now has a variety of interchangeable elbows, which creates multiple uses for just one mask. Along with interchangeable elbows, a selection of cushion shapes are provided to help find the right fit for patients.



Classic
(PD)™



Chin Support
(FP)™



Wide Profile
(RM)™



High Profile
(RP)™



Standard



Anti-Asphyxia



Anti-Asphyxia
with Leak



Standard
Nebulizer



Anti-Asphyxia
Nebulizer



Standard
Bronchoscopy




Anti-Asphyxia
Bronchoscopy



For more information about BiTrac Select™ and our complete line of ventilation products visit our website at www.pulmonary.com or call 317.246.5505.

.... bringing change to life[®]



Part No.	Size	Case/Qty.
Masks (All masks include Head Strap)		
BiTrac Select™ NIV Full Face Mask with Interchangeable Anti-Asphyxia (EE) Elbow 22mm Female and OmniClip™ with Silicone Forhead Pad		
Classic Cushion (PD)™		
313-9504	Adult Extra Large	10
313-9503	Adult Large	10
313-9502	Adult Medium	10
313-9501	Adult Small	10
Chin Support Cushion (FP)™		
313-9563	Adult Large	10
313-9562	Adult Medium	10
313-9561	Adult Small	10
Wide Profile Cushion (RM)™		
313-9533	Adult Large	10
313-9532	Adult Medium	10
313-9531	Adult Small	10
313-9530	Adult Extra Small	10
High Profile Cushion (RP)™		
313-9543	Adult Large	10
313-9542	Adult Medium	10
313-9541	Adult Small	10
BiTrac Select™ NIV Full Face Mask with Interchangeable Standard (SE) Elbow 22mm Female and OmniClip™ with Silicone Forhead Pad		
Classic Cushion (PD)™		
313-9509	Adult Extra Large	10
313-9508	Adult Large	10
313-9507	Adult Medium	10
313-9506	Adult Small	10
Chin Support Cushion (FP)™		
313-9568	Adult Large	10
313-9567	Adult Medium	10
313-9566	Adult Small	10
Wide Profile Cushion (RM)™		
313-9538	Adult Large	10
313-9537	Adult Medium	10
313-9536	Adult Small	10
313-9535	Adult Extra Small	10
High Profile Cushion (RP)™		
313-9548	Adult Large	10
313-9547	Adult Medium	10
313-9546	Adult Small	10
BiTrac Select™ NIV Nasal Mask with Standard Elbow 22mm Female and Silicone Forhead Pad		
Classic Cushion (PD)™		
313-9016	Adult Extra Large	10
313-9015	Adult Large	10
313-9014	Adult Medium	10
313-9013	Adult Small	10
Wide Profile Cushion (RM)™		
313-9039	Adult Large	10
313-9028	Adult Medium	10
313-9030	Adult Small	10
High Profile Cushion (RP)™		
313-9049	Adult Large	10
313-9041	Adult Medium	10
313-9040	Adult Small	10
Interchangeable Elbow Options		
313-9599	Standard (SE) Elbow - for dual limb circuit	10
313-9598	Anti-Asphyxia (EE) Elbow - for single limb circuit	10
313-9597	Standard (SE) Bronchoscopy Elbow - for dual limb circuit*	10
313-9596	Anti-Asphyxia (EE) Bronchoscopy Elbow - for single limb circuit*	10
313-9595	Standard (SE) Nebulizer Elbow - for dual limb circuit*	10
313-9594	Anti-Asphyxia (EE) Nebulizer Elbow - for single limb circuit*	10
313-9593	Anti-Asphyxia With Leak Elbow	10
*Available Soon		
<div> <div>Made in USA</div> <div>Single Patient Use</div> <div>Rx ONLY</div> <div> Latex Free </div> </div>		

