



Trans-vaginal

REF 1204 series



Description

Trans-vaginal and trans-rectal transducers are widely used in gynaecological and urological procedures. This technique brings the transducer into closer proximity to the target organs providing optimal images for use in diagnosis.



Packaging information

Packaging Single packed
Selling unit 05x50 (non sterile)
05x24 (sterile)



Specifications

Sterile	Yes
Non-sterile	Yes
Material	Polyurethane
Latex	No
Powder	No
PVC	No
Phthalates	No



Dimensions

REF 1204-01 30 x 300 mm including twist-lock non-sterile
REF 1204-02 30 x 300 mm including twist-lock sterile

Probetection™ trans-vaginal covers

The Probetection™ ultrasound transducer covers offer an optimal fit for all trans-vaginal transducers currently in use.

All transducer covers are leak tested, single packed and available sterile and non-sterile. The Probetection™ polyurethane transducer covers are successfully tested to be impenetrable for viruses and bacteria.

Probetection™ transducer covers can be easily applied aseptically.

Probetection™ protects

Probetection™ offers the perfect protection for vaginal and obstetric ultrasound transducers. Palmedic Probetection™ ultrasound transducer covers protect the patient, the healthcare worker and the equipment against cross contamination.

Probetection™ helps fighting nosocomial infections, including the transmission of Human Papillomavirus, Hepatitis and other kinds of pathogenic microorganisms.

Each cover comes with a twist-lock to secure it to the cable of the transducer.

